IN THE UNITED STATES PATENT & TRADEMARK OFFICE

Applicant(s): Tereba et al. Docket No.: 016026-9043-US01

Serial No.: 10/694,475 Group Art Unit: 1639

Filing Date: October 27, 2003 Examiner: Christopher M. Gross

Title: SIMULTANEOUS ISOLATION AND QUANTITATION OF DNA

DECLARATION OF REX BITNER UNDER 37 CFR § 1.132

Assistant Commissioner for Patents Washington, DC 20231

Sir:

I, Rex Bitner, do hereby declare and state the following:

- I have served as the Technology Manager of the Genetic Analysis R&D group at Promega Corporation, Madison, WI, since 2003. From 1997-2003, I worked as a Senior Scientist and Senior Project Manager of R&D at Promega Corporation. I hold a B.S. in Biology from The University of Washington, Seattle, WA and a Ph.D. in Genetics from The University of California, Davis, CA. A copy of my curriculum vitae is attached as Exhibit A.
- I am a joint inventor of at least some of the claimed subject matter of the above-identified patent application. I make this declaration in support of prosecution of the present application before the U.S. Patent and Trademark Office.
- 3. I have read and understand the invention as disclosed in the present application, including the invention described by the presently pending claims. I have also reviewed the July 31, 2007 Office Action. I understand that each of claims 44-52, 54, 55, and 58-68 is rejected as being unpatentable (i.e., obvious) over Kleiber et al. (WO 96/41811), Huber et al. (1993 Nuc. Acids Res. 21:1061-1066), and Vogelstein et al. (1979 PNAS 76:615-619).
- Claims 44-52, 54, 55, and 58-68 are directed to methods of isolating a defined and consistent amount of DNA from multiple samples by choosing the amount of DNA to be

isolated, choosing an amount of a silica containing solid support needed to isolate the defined amount of DNA, such that the amount of DNA in the samples is greater than the binding capacity of the solid support, and contacting each sample with the solid support under conditions that allow isolation of the defined and consistent amount of DNA.

- 5. DNA IQTM System is the term that Promega Corporation, the assignee of the instant application, uses with its customers when referring to methods for isolating a defined and consistent amount of DNA from multiple samples, as described throughout the application and as summarized in paragraph 4, above.
- 6. The concept underlying the methods of the invention, which focus on isolating just a portion of DNA that may be present in a sample, represents a complete departure from prior art methods such as those described in Kleiber et al. and Vogelstein et al., which focus on maximizing DNA yield. For example, Kleiber et al. discusses the relatively high yields obtained by their methods. (Please see, for example, page 12 of Kleiber et al.). Similarly, Vogelstein et al. emphasizes that binding DNA to glass from dissolved agarose "is rapid, convenient, and nearly quantitative." (Vogelstein et al., p. 618, first column, last line to second column, first line, emphasis added). Huber et al., rather than isolating a defined and consistent amount of DNA, describes using high performance liquid chromatography (HPLC) to fractionate relatively small (i.e., <500 bp) DNA fragments such as restriction fragments and PCR products according to size.</p>
- 7. The ability to isolate a defined and consistent amount of DNA from multiple samples using the methods of the invention simplifies processing, reduces the amount of time needed to process samples, and increases sample throughput. Because DNA samples prepared using the method of the invention contain a select defined and consistent amount of DNA, the isolated DNA can be used directly in downstream applications requiring an amount of DNA within a particular range.
- 8. The importance of increased sample throughput in isolating DNA cannot be overstated. For example, because prior art methods of isolating DNA from samples were very labor-intensive and cost about \$1000 per test, a backlog of several hundred thousand samples from rape victims awaited processing in 1999 when the patent application for the present invention

was submitted. Thus, samples containing evidence that may have been useful in identifying sex offenders go unprocessed.

- Attached as Exhibit B is an article entitled "Robotic Extraction of Mock Sexual Assault Samples using Biomek® 2000 and the DNA-IQTM System", authored by Susan Greenspoon and Jeff Ban of the Virginia Division of Forensic Science, which appeared in the February 2000 issue of Profiles in DNA.
- 10. Greenspoon and Ban describe how the DNA IQTM System, used in conjunction with the Biomek® 2000 robotics system, was able to isolate uniform amounts of DNA from mock sexual assault samples containing dilutions of semen of from 1:10 to 1:200 on a ½, ¼, or 1/8 portion of a swab (See p. 4, Mock Sexual Assault Samples). The DNA thus isolated from samples having widely varying DNA content was used directly in PowerPlex® 1.1 System, a DNA amplification reaction used in genetic identity testing, and produced uniform results.
- 11. Other experts in the field have recognized the importance of the instantly claimed methods isolating DNA for use in molecular biological methods such as amplification for genetic identity testing.
- 12. Dr. Weimin Sun, Scientific Director in the Molecular Genetics Department of Quest Diagnostics Nichols Institute, San Juan Capistrano, CA, evaluated the DNA IQTM System for use in clinical samples and found that the methods yielded consistent quantities of DNA despite significant variations in the starting material, which afforded satisfactory performance in downstream applications. (See correspondence form Dr. Sun to Mr. David Phelps, Exhibit C).
- 13. Kim Gorman, then President of Paternity Testing Corporation, Columbia, MO, reported that the DNA IQTM System was used to extract DNA from buccal swabs in a 96 well format. (See correspondence from Ms. Gorman to Mr. Phelps, Exhibit D). Ms. Gorman noted that buccal swabs vary greatly in DNA content, and also noted that DNA multiplexes (amplification of multiple DNA loci in a single reaction) are concentration sensitive. Despite these challenges, Ms. Gorman further noted that there was no need to quantify the DNA prepared using the DNA IQ extraction method prior to use in a multiplex reaction. Ms. Gorman reported that, using the DNA IQTM System, the total time required to process 96 samples was reduced from 6 or 7 hours

using their traditional extraction method to about 3 hours, and hands on time by analysts was reduced from more than 4 hours to less than 30 minutes.

- 14. Jeffrey Ban, Section Chief of Forensic Biology at the Virginia Department of Criminal Justice Services, Division of Forensic Science, reported that, using the DNA IQTM System in conjunction with Beckman Coulter Biomek® 2000 Workstation or other similar instrumentation, could greatly increase a laboratory's throughput capabilities, thus permitting the Forensic DNA community to provide better service to law enforcement agencies (See correspondence from Mr. Ban to Mr. Phelps, Exhibit E).
- 15. In addition, Promega has been contacted by numerous law enforcement agencies that ultimately used the DNA IQTM System to analyze forensic samples from crime scenes. For example, the Royal Canadian Mounted Police requested assistance in processing samples collected from a hog farm in British Columbia, where the partial remains of at least 26 murder victims were found. The DNA IQTM System was also used to isolate DNA for genetic identity testing of victims at Ground Zero and human remains found in mass graves in Bosnia.
- 16. In 2002, Promega received an R&D 100 Award for its DNA IQ[™] System. Through the R&D 100 Awards program, sponsored by R&D Magazine, organizations receiving the awards are recognized for the most technologically significant products introduced into the marketplace.
- 17. I hereby declare that all statements made herein of my own knowledge are true and that all statements made on information and belief are believed to be true; and further that these statements and the like so made are punishable by fine or imprisonment, or both, under Section 101 of Title 18 of the United States Code, and that such willful false statements may jeopardize the validity of the application or any patent issuing thereon.

Dated: October 30, 2007	RexBitner
	Rex Bitner

O:\CLIENT\016026\9043\B1171888.1

REX M. BITNER, Ph.D.

W53 N598 Birch Street Cedarburg, WI 53012 Work: 608-277-2491 rex.bitner@promega.com

QUALIFICATIONS SUMMARY

Over twenty years of industrial experience in biotechnology research and product development at Promega (Madison, WI), Amersham/Pharmacia Biotech (now GE Healthcare), (Milwaukee, WI), and 3M (St. Paul, MN). Four years of postdoctoral research in molecular biology at the University of Colorado, Boulder and the University of California, Davis. Broad, in-depth knowledge of molecular biology, solid phase purification of biological materials, and all aspects of product development in an ISO9001: 2000 environment. Experienced in technology management, project leadership and management, supervision of a Bl.3 laboratory, and the development, launch and care for biotechnology products. Extensive experience in the creation, development and evaluation of intellectual property, including USPTO and foreign patent office actions, in person interviews with USPTO patent examiners, as well as participation in foreign patent opposition proceedings, both as proprietor and as an opposition party, and as an opposition member in the appeal of an invalidated patent.

PROFESSIONAL HISTORY

1997 - present PROMEGA CORPORATION, Madison, Wisconsin

Technology Manager, Genetic Analysis, R&D (2003 - present)

Technology Manager, Genetic Analysis: management of a laboratory group, development of separation technologies, intellectual property, and purification products for the biotechnology marketplace and clinical laboratory market, with particular emphasis on automation of nucleic acid purification products for use in genomics, high throughput pharmaceutical drug screening, and clinical diagnostics. Intellectual property management has included USPTO and foreign patent office actions, in person interviews with USPTO patent examiners, phone interviews with USPTO and foreign patent examiners, and participation in European patent opposition proceedings, both as proprietor and as an opposition party, including as an opposition member in the appeal of an invalidated patent. The efficient integration of foreign patent filling strategies and licensing of intellectual property within Promega's Genetic Analysis business strategy has been a central responsibility.

Senior Scientist, Senior Project Manager, R&D (1997 - 2003)

Development of pH dependent ion exchange purification systems using both column and paramagnetic particle purification methods in robotic workstations, with additional emphasis on automated cell concentration and magnetic clearing of cellular lysates in 96-well walkaway automated DNA purification (particularly using Beckman BioMek® FX and Tecan Genesis® robotic platforms). Products developed for genomic DNA purification from human whole blood and tissues (including R&D Magazine's R&D 100 award winning DNA-IQ (for 2002), plant materials, and DNA purification from food ingredients for use in the quantitative detection of genetically modified organisms (GMO) in food. Additional experience with DNA sequencing automation, RNA purification, and PCR cleanup, particularly using Beckman, Tecan and Thermo-Electron Labsystems robotic platforms.

Project leader for Promega products, in an ISO9001: 2000 environment, including:

Wizard® Genomic, 10ml blood

A1620

Wizard® Magnetic DNA Purification for Food FF3751

Wizard® Magnesii® Dlood Genomic, Max Yield

MagneSii® NOE, Fixed Yield

MD1370

MagneSii® KF, Genomic System

MD1460

PureYield™ Plasmid Midi-Prep System

A2495

PROFESSIONAL HISTORY, continued

1995 - 1997 AMERSHAM PHARMACIA BIOTECH INC., Milwaukee, Wisconsin.

Senior Research Scientist (1995 - 1997)

Senior Research Scientist and Project Leader responsible for nucleic acid purification products for the biotechnology laboratory: development of novel separations matrices and processes, formulation and execution of project plans, maintenance of timelines and scheduling, and supervision of personnel on several project teams. Research and development of new products for the molecular biology marketplace, particularly in the areas of proprietary purification products, anion exchange chromatography, and solid phase extraction and immobilization of nucleic acids. Responsibilities included the management of personnel, timelines and ISO 9001 documentation of product development.

1982 - 1994 3M COMPANY, St. Paul, Minnesota.

Research Specialist (1984 - 1994)

Identifying, planning and pursuing molecular biology programs of interest to 3M businesses. Research programs involved diverse product objectives: Development of 3M's Rapid AttestTM biological/Sterilization monitor product, GMP purification of bovine phosphophoryn proteins (for bone repair), genetic manipulation of bacteria to produce specialty chemicals (meta-hydroxyphenylacetylene, and aromatic compounds useful in laser dyes), cDNA cloning of mammalian genes for drug discovery screening, surface immobilization of nucleic acids onto ceramic oxide/3M EmporeTM (PTFE) membranes (for use in DNA blotting, hybridization and sequencing), solid phase extraction of DNA for automated sequencing, solid phase extraction of DNA from human blood plasma for use in PCR, protein immobilization on azlactone functionalized porous beads (including 3M EmphazeTM beads and network beads), and cloning of stress protein genes from bacteria associated with periodontal disease.

Development and implementation of DNA purification technologies, using a variety of ceramic matrices for solid phase extractions. Over ten years experience supervising a biosafety level 3 containment laboratory.

Other responsibilities included evaluation of both internal and external research proposals and intellectual property. Additional responsibilities as Institutional Biosafety Officer, OSHA blood-borne pathogen safety officer, and as scientific advisor in the development of and compliance with Minnesota State (Environmental Quality Board) regulations governing recombinant organisms.

Senior Biologist (1982 -1984)

Responsibilities included setting up and staffing a recombinant DNA laboratory, evaluation of outside business proposals and intellectual property issues, and initiation of new research programs in molecular biology: gene expression in Bacillus subtilis, and R&D of 3M's Rapid Attest. Additional responsibilities: initiation of university research contracts and management of laboratory personnel.

1978 - 1982 UNIVERSITY OF COLORADO, Boulder, Colorado.

Postdoctoral research: Dr. Peter L. Kuempel, Dept. of Molecular, Cellular, and Developmental Biology: Termination of chromosome replication in *E. coli*.

PROFESSIONAL HISTORY, continued

1974 - 1978 UNIVERSITY OF CALIFORNIA, Davis, California.

Postdoctoral research: Dr. Gordon G. Edlin, Dept. of Genetics (1978). Instructor: Department of Genetics (1977).

1974 UNIVERSITY OF WASHINGTON, Seattle, Washington Post-graduate Research Assistant: Dr. Jonathan A. Gallant. Dept. of Genetics

EDUCATION

Ph.D. in Genetics, 1978 University of California, Davis, California

B.S. in Biology, cum laude, 1974 University of Washington, Seattle, Washington

PATENTS

Rex M. Bitner, Mark A. Denhart, Don Smith. US Utility Application: WO2007/103485. Small RNA Purification.

Rex M. Bitner, Paula R. Brisco, Michelle N. Mandrekar. WO2007070381, US20050748825. Nucleic Acid Purification with a Binding Matrix.

Rex M. Bitner, Michelle Mandrekar, Don Smith, Douglas H. White. WO2007005613, US2007015191.
Network of Buoyant Particles for Biomolecule Purification and Use of Buoyant Particles or Network of Buoyant Particles for Biomolecule Purification.

Rex M. Bitner, Daniel J. Simpson, Roderick G. Flemming and Susan C. Koller. US 6,787,307, EU 1341910, CA 2428532, AU 2594202 and WO0238758. Lysate Clearance and Nucleic Acid Isolation using Silanized Silica Matrices.

Allan M. Tereba, Rex M. Bitner, Susan C. Koller, Craig E. Smith, Daniel D. Kephart, Steven J. Ekenberg. US 6,673,631 B1, CA 2379503, EP 1204741, EP 1510578, DE 60015148 and WO0114590. Simultaneous Isolation and Ouantitation of DNA.

Rex M. Bitner, Jacqui Sankbeil, Braeden L. Butler, Douglas H. White, Craig E. Smith. US 7,078,224, US 6,284,470, EP1179058, EP1341910, EP1621618, and WO0070040. Cell Concentration and Lysate Clearance Using Paramagnetic Particles.

Craig E. Smith, Diana L. Holmes, Daniel J. Simpson, Jehoshua Katzenhendler, Rex M. Bitner, Josephine C. Grosch. US 6,806,362, US 6,310,199, EP1179057, AU5126100 and WO0069872. pH Dependent Ion Exchange Matrix and Its Use in the Isolation of Nucleic Acids.

Craig E. Smith, Diana L. Holmes, Daniel J. Simpson, Jehoshua Katzenhendler, Rex M. Bitner, Josephine C. Grosch. US 6,376,194, US 6,270,970, EP1179056, AU4841500 and WO0070041. Mixed Bed Solid Phase and Its Use in the Isolation of Nucleic Acids.

PATENTS, continued

Craig E. Smith, Donald A. Creswell, Rex M. Bitner, Douglas H. White, Braeden L. Butler, Scott A. Lesley. US 6,194,562, EP1071695, CA 2329067, and WO 9954340. Endotoxin Reduction in Nucleic Acid Purification.

Rex M. Bitner, Chan-Wha Kim, and Michael G. Williams. EP647232B1, WO9400464. 3M, St. Paul, MN. Deproteinization with azlactone-functional supports.

Rex M. Bitner and Eric F. Funkenbusch. EP0391608. Applicant: Minnesota Mining and Manufacturing Company, St. Paul, MN. Metal oxide supports for nucleic acids.

PUBLICATIONS

Koller, S., H. Shenoi, R. Bitner, 2005. "Purifying Genomic DNA from Plant Tissue on Automated High-Throughput and Moderate-Throughput Platforms" JALA 10: 155-162.

Bitner, R. M. and S.C. Koller. 2004. "Automated high throughput purification of genomic DNA from plant leaf or seed using Promega's MagneSil® paramagnetic particles? "Proceedings of SPIE, Microarrays and Combinatorial Techniques: Design. Fabrication and Analysis II. Vol 5328. p. 78-86.

Bitner, R., S. Koller, J. Sankbeil, M. Denhart and H. Shenoi. 2004. "Purifying Genomic DNA from Whole Blood on Automated, High-throughput and Moderate-throughput Platforms" JALA 9: p. 64-71.

Smith, C., P. Otto, R. Bitner and G. Shiels. 2003. "Chapter 9: DNA Purification", in <u>PCR Primer: A Laboratory Manual</u> 2nd Edition, editors: Carl W. Dieffenbach and Gabriela S. Dveksler, Cold Spring Harbor Laboratory Press, ISBN 0-87969-6532, p.87-115.

Bitner, R., S. Koller, and J. Sankbeil. 2003. "Automated high throughput purification of genomic DNA from whole blood using Promega's MagneSil® paramagnetic particles with either the Max Yield or MagneSil® ONE normalized purification methods "in Microarrays and Combinatorial Technologies for Biomedical Applications: Design, Fabrications and Analysis, Proceedings SPIE 4966, p. 98-105.

Bitner, R.M. and S.C. Koller. 2002. "Automated genomic DNA purification options in agricultural applications using MagneSil" paramagnetic particles. Tools for Molecular Analysis and High-Throughput Screening Proceedings of SPIE Vol. 4626 p. 218-225.

Bitner, R. and S. Koller. 2001. "Automation of DNA Extraction from Food and Plants Using MagneSilTM Paramagnetic Particles" Genomics and Proteomics Technologies. *Proceedings of SPIE* Vol 4264 p. 9-16.

Bitner, R., D. White, S. Krueger, M. Bjerke, B. Butler, C. Smith. 2000. "Use of MagneSil™ Paramagnetic Particles for Plasmid Purification, PCR Cleanup and Purification of Dideoxy and Big Dye DNA Sequencing Reactions" Advances in Nucleic Acid and Protein Analyses, Manipulation and Sequencing, Proceedings of SPIE Vol 3926 p. 126-133.

Williams, M.G., P.E. Olson, K.J. Tautvydas, R.M. Bitner, R.A. Mader, and L.P. Wackett. 1990. "The application of toluene dioxygenase in the synthesis of acetylene terminated resins" Appl Microbiol Biotechnology 34: 316-321.

Bitner, R.M. and P.L. Kuempel. 1982. "PI Transduction Mapping of the trg Locus in rac* and rac Strains of Escherichia coli K-12" J. Bacteriol. 149: 529-533.

REX M. BITNER, page 5

PUBLICATIONS, continued

Bitner, R.M. and P.L. Kuempel. 1981. "P1 Transduction Map Spanning the Replication Terminus of Escherichia coli K-12" Mol. gen. genet. 184: 208-212.

Binding, R., G. Romansky, R. Bitner and P. Kuempel. 1981. "Isolation and Properties of Tn10 Insertions in the rac Locus of Escherichia coli" Mol. gen. genet. 183: 333-340.

Edlin, G., L. Lin and R. Bitner. 1977. "Reproductive Fitness of P1, P2 and Mu Lysogens of Escherichia coli" J. Virol. 21: 560-564.

Lin, L., R. Bitner and G. Edlin. 1977. "Increased Reproductive Fitness of Escherichia coli Lambda Lysogens" J. Virol. 21: 554-559.

Gallant, J.A., L. Shell and R. Bitner. 1976. "A Novel Nucleotide Implicated in the Response of E. coli to Energy Source Downshift" Cell 7: 75-84.

PROFILES IN DNA

Volume 5, No. 1

February 2002

In This Issue

Automated Systems for DNA Analysis

DNA Database Legislation and Legal Issues

Development of PowerPlex® Matrix and Sample Protocols on the ABI PRISM® 3100

DNA IQ™ System FAQs

Robotic Extraction of Mock Sexual Assault Samples Using the Biomek® 2000 and the DNA IQ[™] System

By Susan Greenspoon and Jeff Ban Virginia Division of Forensic Science, Richmond, Virginia

INTRODUCTION

Forensic scientists are routinely faced with the challenge of isolating DNA from a large array of tissue and cell types. The variety of substrates upon which colladar material has been deposited, some of which may contain inhibitous of PCR, can make the process more difficult (1). Therefore, any sobotic system applied to the extraction of forensic casework sampler must be robust enough to address these variations. The Biomét '2000 used in conjunction with the DNA IQP' System⁶⁰ has proven to be an efficient robotic system designed to handle the challenges of routine casework samples.

The DNA IQ²⁰⁰ System uses silica-coated magnetic beads to separate DNA from cellular debrits. Cells are lysed in a powerful lysis buffer, and the lysate is mixed with the magnetic beads. The beads saturate at approximately 100ng of bound DNA, and the excess DNA is removed by pipetting. Once bound to the magnetic resin, the DNA is pipetted and vigorously shaken several times in wash buffer, then eluted using heat. The Biomel* 2000 is equipped with a magnetic plate, a shaking platform and a thermal exchange unit to perform these necessary steps.

CONTAMINATION STUDIES

A number of exploratory and validation studies have been performed on the Biomek* 2000.
DNA IO* System to evaluate the viability of this automated system for true with forenzie samples "Manamental questions needed to be answered before moving ahead with extensive validation work. First, does the robotic, open-plate format cause contamination? To narver this question, two sample formats were used repeatedly, and extracted amplies were analyzed to test

The Biomek® 2000 used in conjunction with the DNA IQ™ System has proven to be an efficient robotic system designed to handle the challenges of routine casework samples.

CSF	1 2 Walled	3 4 MIII	1 5 6 1 1 2	7 6 =	9	10 11 =	12 13 ======	14 +• """"		_	_				=	-	
трох	1111			F		1	-111	1111111	D7	=	Ē	- - - -	1111111	=	milin	-	
TH01	151111			-	111111	-1	=	1	D13	=======================================		-	111111111	4		The state of the s	
vWA	iiillinu i	ndanih i	milding i					i indilan	D			million		111111111111111111111111111111111111111			METATO.IA

Figure⁴. A Amplified DNA samples using the PowerPlee⁴. 1.1 System from the 88-sample checkerboard contamination study. Numbers 1-14 indicates amplier w² indicates positive control (KGS)3, and ²-indicates applier control. Amplified DNA was analyzed by decruptonesis using a 6% polyacrylamids gd (Ghco BRL) for 2 hours at 50 wats. Both the SSS mis scan (left panel) and the 505 mm scan (right panel) are shown (gd imaging performed using a Blackin Bladfor instrument).

Profiles in DNA/February 2002

FEATURE ARTICLE

for contamination. The first is the zebrastripe format test; alternating columns of samples containing an abundant source of DNA with columns containing reagent blanks (8 sample wells per column). Therefore, a column of samples containing abundant DNA was processed adjacent to a column of reagent blanks in a striped pattern on the plate. The samples containing DNA were bloodstains cut into 5mm2 squares. DNA was eluted from the magnetic beads into 100ul of sterile water, then quantified, amplified, and typed using the PowerPlex 1.1 System(be) (2). The first two trials of this test detected some contamination. The software method used was modified to accommodate sample loading into a 96 deepwell plate in place of the more shallow Greiner plate and to remove an initial shaking step. A subsequent zebra stripe experiment showed no contamination with the 40 samples that were isolated.

The second contamination test is a checkerboard sample format: samples containing abundant DNA were alternated with reagent blanks in a checkerboard pattern across a 96 deep-well plate (Figure 1). All 128 samples (88 sample and 40 sample methods) tested negative for any detectable contamination.




Figure 2. PowerPlex" 1.1 mock sexual assault comparison study. Duplicate samples were extracted manually (indicated by an "O" above the lane, for organic) or robotically (indicated by a "B" above the lane, for Biomek* 2000). Semen dilutions (1:10, 1:100 and 1:200) for each set of swabs are indicated above the corresponding six sample lanes, which contain DNA extracted from the indicated swab portion (1/2, 1/4 or 1/8). Amplified DNA was separated by electrophoresis in a 6% polyacrylamide gel (Gibco BRL) for 2 hours at 50 watts. Both the 585nm scan (left panel) and the 505nm scan (right panel) are shown (gel imaging performed using a Hitachi FMBIO* instrument).

MOCK SEXUAL ASSAULT SAMPLES

Sexual assault cases frequently constitute the majority of DNA cases received by a forensic laboratory, Presently, no robotic system is available that can separate sperm from non-sperm cells and thus perform a differential extraction (3) from start to finish. However, the first step of separating fractions can be performed manually. Subsequently, the DNA from E-cell (epithelial or nonsperm cell) lysates and sperm pellets can be extracted robotically, saving analysts a substantial amount of time. Any robotic extraction of sexual assault samples must at least be able to generate sample DNA, equivalent in both quality and yield to that generated by manual DNA extraction methods. Therefore, a thorough examination of the Biomek® 2000/DNA IQ™ System's ability to isolate DNA from sexual assault samples needed to be performed. The first step was to ascertain whether sperm cells could be successfully lysed and the DNA purified by the robotic system. Mock sexual assault samples were prepared using previously donated vaginal swabs and semen from a known donor, which was deposited onto sterile cotton swabs in 1:2 and 1:4 dilutions. The E-cells were lysed manually, the sperm cells pelleted, and a portion of the

lysates and the entire sperm pellets were loaded onto the Biomek" 2000 for DNA extraction. DNA was eluted off the magnetic beads into 100µl of sterile water. Highquality DNA was obtained and typed accurately using the PowerPlex* 1.1 System (data not shown).

Once it was demonstrated that the Biomek® 2000/DNA IO™ System could successfully complete the differential extraction process with the E-cell lysates and intact sperm, the next question addressed was whether this automated system could produce DNA of comparable quality and yield to that produced by manual extraction. A comparative study was designed to measure the performance of the Biomek" 2000/DNA IQ™ System with respect to manual extraction of similar if not identical samples. Samples were prepared in the following manner:

- 1. Sets of vaginal swabs from 5 different donors were selected.
- 2. Duplicate mock sexual assault swabs were prepared using semen from a single donor at the following dilutions: 1:10, 1:100, 1:1,000 and 1:10,000 for three sets, 1:10, 1:100, 1:200 and 1:400 for one set and 1:100, 1:200, 1:400 and 1:800 for the
- 3. Once dried, the swabs were cut into 1/2, 1/4 and 1/8 portions.
- 4. The E-cells were lysed and the sperm cells pelleted and washed.
- 5. The samples were split evenly with one half going to an analyst to complete the extraction manually and the other loaded onto the Biomek* 2000 for a robotic DNA extraction

Yields and quality of the DNA from the sperm fractions processed by the Biomek® 2000/DNA IQ™ System were comparable and frequently superior to those obtained by manual extraction (Figure 2). In less experienced hands, the Biomek® 2000/DNA IO™ System clearly outperformed the manual extraction (data not shown). Results obtained by experienced users were equivalent to those achieved with the robot. Therefore, the Biomek* 2000/DNA IO™ System is not only capable of outperforming its human counterpart, but it also delivers a more consistent product.

www.promega.com

The maximum sample volume for use with the 96 deep-well plate is limited to 100µL Sperm cells are typically in a pellet of 50µl and therefore unaffected by the volume limit. Since the E-cell lysate is usually in a volume of 500µl, it was important that the yields from 100µl of E-cell lysate from the 1/2, 1/4 and 1/8 swab portions be sufficient for all DNA typing needs. Total yields of E-cell DNA extracted on the robot were calculated (Figure 3), and sufficient E-cell DNA could be obtained using robotic extraction methods. In fact, 1/8 swab portions provided more than enough DNA for an E-cell DNA profile.

The amount of time saved by using the Biomek* 2000/DNA IQ™ System to extract the forensic samples can be substantial. The time it takes to complete the organic extraction manually for a single sample is 5 hours and 5 minutes (after E-cells have been lysed, and sperm cells pelleted and washed). Of course, additional samples will lengthen the amount of time proportionately. In comparison, the robot takes 1 hour and 15 minutes to extract the DNA from 40 samples and 1 hour and 50 minutes to extract 88 samples. Therefore, the absolute minimum amount of time saved is 3 hours and 50 minutes or half of a day.

EXTRACTION OF OTHER CELL AND TISSUE TYPES

Since a variety of cells and tissue types are encountered in routine forensic casework, the Biomek* 2000/DNA IQ™ System was evaluated to determine its capability to isolate DNA from a variety of sources. Dried bloodstains, E-cell lysates, intact sperm cells, muscle, heart, brain, liver and buccal swabs were extracted using the Biomek® 2000/DNA IQ™ System and successfully typed using the PowerPlex* 1.1 System (data not shown).




Figure 3. Bar graph depicting the total yield of E-cell DNA generated from extraction on the Biomek* 2000

robot using 100µl of lysate. Yields were determined by measuring DNA concentration using the QuantiBlot kit.

CONCLUSION

A completely automated system for extraction of sexual assault samples is currently not available. However, once the E-cells have been separated from the sperm cells, robotic DNA extraction using the Biomek* 2000/DNA IQ™ System can be accomplished. Moreover, when sperm DNA is limited, the Biomek® 2000/DNA IQ™ System generates DNA of similar, and sometimes better, quality and yield than that obtained by manual extraction of a duplicate sample. Because of the adaptability of the Biomek* 2000/ DNA IQ™ System, this instrument has the potential to handle future applications of emerging cell separation technologies. It may be possible on a robotic platform, to separate sperm cells from non-sperm cells with the use of an anti-sperm antibody conjugated to a magnetic bead. One can envision a completely automated system where both cell separation and DNA extraction are performed on the same robot. The Biomek* 2000/DNA IQ™ System may be uniquely poised to proceed with that application when the technology becomes available.

The time saved when compared with manual extraction, as well as the ability to extract a variety of tissue and cell types, makes the Biomek® 2000/DNA IQM System attractive for casework applications. Further validation work on the Biomek® 2000/DNA IQ™ System must be performed in order to complete our evaluation and validation prior to its application for forensic casework. Although no contamination of the samples was detected after modifying the method and changing the format, special circumstances might require the use of manual extraction methods, for example, when an evidentiary sample may be completely consumed due to limited available material

ACKNOWLEDGMENTS

We would like to acknowledge the participation of other members of the Virginia Division of Forensic Science: Beth Ballard, Missy Baisden, Brian Covington, Shelley Smith and Colleen Young, We would also like to thank Allan Tereba and Dan Kephart at Promega Corporation for all their hard work at making the technology and this study possible.

REFERENCES

- 1. Cattaneo, C. et al. (1997) Comparison of three DNA extraction methods on bone and blood stains up to 43 years old and amplification of three different gene sequences. J. Forensic Sci. 42, 1126.
- 2. Lins, A. et al. (1998) Development and population study of an eight-locus short tandem repeat (STR) multiplex system. J Forensic Sci 43, 1.
- 3. Gill, P., Jeffreys, A.J. and Werrett, D.J. (1985) Forensic application of DNA 'fingerprints'. Nature 318, 577.

PowerPlex is a trademark of Promega Corporation and is registered with the U.S. Patent and Trademark Office. DNA IQ is a trademark of Promega Corporation.

Biomek is a registered trademark of Beckman Instruments, Inc. FMBIO is a registered trademark of Hitachi Software Engineering Company, Ltd. QuantiBlot is a registered trademark of Roche Molecular Systems, Inc.

(a,b,c)Refer to the patent and disclaimer statements on page 2.



March 4, 2002

David Phelps Promega Genetic Identity 2800 Woods Hollow Drive Madison, WI 53711

Dear Mr. Phelps,

We have evaluated DNA IQ system for isolation of DNA from clinical specimens (whole blood and buccal swabs). The IQ system showed satisfactory performance on both types of samples. It is highly tolerant to the variability of sample input quantity, It also demonstrated high consistency in the quantity of product DNA. This feature is extremely desirable to its as buccal swab samples can have significant variations in starting material affected by collection procedure and skills. In addition, DNA samples isolated using the DNA IQ system have shown satisfactory performance in the downstream applications and have excellent stability as well.

If there is any question in regards to our experience with the DNA IQ system, please do not hesitate to contact me.

Sincerely

Weimin Sun, Ph.D., ABMG. Scientific Director Molecular Genetics Department Quest Diagnostics Nichols Institute 33608 Ortega Highway

San Juan Capistrano, CA 92690 (949) 728-4498 (voice) (949) 728-4874 (fax)



Paternity Testing Corporation Fast, Confidential DNA Analysis

February 25, 2002

David Phelps Promega, Genetic Identity 2800 Woods Hollow Drive Madison, WI 53711

VIA FACSIMILE: (608) 273-6455

Dear Mr., Phelps:

Lisa Lane requested that I send you a brief letter about our experience using the DNA IQ extraction method.

PTC has been extremely interested in automating DNA extractions from buccal swabs. Because of the nature of buccal swabs it has historically been extremely difficult if not impossible to work with them in a 96 well automated format. Swabs shavehed liquid and had to be spun through a basket in order to recover the lysis buffer. Swabs had to be removed by hand one at a time. There were many tamsfers in and out of centrifuges. It was a very labor intensive process to carry out the extractions. It was possible to automate the extractions, but that required large volumes and could not be carried out in a 96 well format.

Our goal was to be able to put a buccal swab into a 96 well tray and never have to touch it again. DNA IQ allows for extractions to be carried out in the 96 well format from the time the specimens are placed into a 96 well deep well tray until the DNA is transferred to the amplification tray.

DNA IQ also has the advantage of climinating the need to quantify the DNA. Buccal swabs vary greatly in DNA content. Since the DNA multiplexes are somewhat concentration sensitive, it was necessary to get a rough quantification of the amount of DNA present then diluting to appropriate volumes. There is no need to quantify the DNA using the DNA IQ extraction method. The concentration of DNA is consistent from well to well. By using 1µl of the isolated DNA we are able to get consistent results from virtually all specimens.

DNA IQ is also much faster than most extraction methods. It takes approximately 3 hours from the time the swabs are placed into the tray until the DNA is pipetted into the amplification tray. There is less than 30 minutes of analysts time involved in the entire process. Using our traditional extraction method it takes 6 or 7 hours to prepare 96 samples for amplification and more than 4 hours of actual hands on time by the analysts.

Our lab has nearly completed the final stages of testing and validation of DNA IQ. We plan to be online using DNA IQ in March.

Yours truly,

Kin Doma

Kim Gorman President



COMMONWEALTH of VIRGINIA

DEPARTMENT OF CRIMINAL JUSTICE SERVICES

DIVISION OF FORENSIC SCIENCE CENTRAL LABORATORY A Nationally Accredited Laboratory

February 25, 2002

P.O. BOX 999 RICHMOND, VIRGINIA 23218 (804) 786-4707

David Phelps Promega Corporation 2800 Woods Hollow Road Madison, WI 53711-5399

Dear Mr. Phelps,

In order to aid the DNA examiner in the extraction of casework samples, the Virginia Division of Forensic Science, Forensic Biology Section, for the past 6 months has worked with scientists and engineers from the Promega Corporation and Beckman Coulter to evaluate, modify and validate the Beckman Coulter Biomeke 2000 Workstation in conjunction with the Promega DNA [QTM Bolation System. Currently each DNA examiner spends between 6 and 8 hours purifying the DNA from an average sexual assault case consisting of 3 to 5 evidence samples [Because of the adaptability of the Promega DNA [QTM Isolation System the DNA purification can be automated permitting the user to walk away once the samples are setup in the robot.

In the time it would takes 4 DNA case examiners (i.e. between 24 and 32 person hours) to manually purify the DNA from evidence samples from 3 to 5 cases, using an organic extraction procedure, the robot can accomplish in under 2 hours due to the speed of the Promega DNA IQT purification process paired with the Beckman Coulter Biomek* 2000 Workstation. Therefore, with as little as 5 to 10 minutes of a technician's time to setup the instrument, the Virginia Division of Forensic Science, Forensic Biology Section can improve the turnaround time of cases with minimum impact on the casework examiner.

Utilizing the chemistry of the Promega DNA IQ™ Isolation System in combination with the Beckman Coulter Biomek® 2000 Workstation or similar instrumentation can greatly increase a laboratory's throughput capabilities for analyzing forensic casework samples. Thus in turn

Fax Numbers:

Administration (804) 786-6857 * Breath Alcohol and Training (804) 786-6139 * Evidence Receiving (804) 786-6924 * Drugs, Trace Evidence Toxicology & Firearms (804) 786-5985 * Latent Prints (804) 786-6190 * Maintenance (804) 786-5985 * Latent Prints (804) 786-6190 *

permitting the Forensic DNA community to provide a better service to the law enforcement agencies throughout the United States to help solve crimes against persons and property.

Yours sincerely,

Jeffrey D. Ban Forensic Biology Section Chief